


KEN-KAGE



Welded one-piece 22 gauge stainless steel body

Stainless steel truss head machine screw

Coved rear corners, top and bottom

Stainless steel, electropolished removable door, can be mounted to open from the left or right side

12 gauge stainless steel latch brackets adjustable for smooth operation

Stainless steel frame

Latch connector notched for easy opening

Stainless steel crossmembers

Stainless steel vertical spokes with 25 mm space between spokes

12 gauge stainless steel self-locking latch connector designed to accommodate a padlock

Stainless steel removable hinge pins

Stainless steel latch keepers

12 gauge stainless steel die formed barrel style hinge

18 gauge stainless steel die formed floor support braces

Door mounts on the front of the cage leaving a minimal gap between the door and cage body for greater animal safety

Door extends below cage floor to retain paper liners

Fluid retaining lip, standard on all cages

KEN-KAGE